
Difficult Airway Management Simulator



The D.A.M. Simulator provides a combination of controllable degree of mouth opening, changeable limit of neck flexibility, and the inflatable tongue.

This model offers wide variations of experiences training with the Difficult Airway Management training.

• 24 variations (includes 1 normal case) of settings by combinations of: 2 steps of neck flexibility, 3 steps of mouth opening, 2 sizes of tongues and 2 positions of vocal codes

Airway Opening Techniques (Head-Tilt/Chin-Lift, Jaw Thrust)


The neck's flexibility can be changed in two steps. The degree of the mouth opening can be set in three steps.





Tongue Edema

The tongue can be inflated to represent the swollen tongue. Laryngospasm Realistic laryngospasm





Laryngospasm

Airway Opening Techniques

True-to-life articulation allows for head-tilt/ chin-lift and jaw-thrust techniques.





Successful ventilation can be confirmed by indicators.






Intraoral Intubation with Laryngoscope

Setting the head at "sniffing position", intubation with laryngoscope can be performed.

Success-Failure Confirmation Successful tube tip placement can be

confirmed by indicators.





Use of Video Laryngoscopes



Use of **Laryngeal Mask**



Training Skills

- -Airway opening techniques (head tilt, jaw thrust)
- -Bag-Valve-Mask ventilation
- -Pre intubation airway assessment
- -"Sniffing position"
- -Pressurization of external larynx to improve the laryngeal view
- -Intraoral/Intranasal Intubation
- -Use of oropharyngeal airway (OPA)
- -Use of nasopharyngeal airway (NPA)

- -Use of laryngeal mask airway
- -Use of video laryngeal scope
- -Use of tracheal intubation fiber scope
- -Confirmation of successful ventilation by indicators.
- -Feedbacks on incorrect procedures including esophagus intubation and unilateral intubation
- -Securing the tube in place with tapes or Thomas™ endotracheal tube holder

Set Includes

- 1 manikin
- 3 upper incisors
- 1 lubricant
- 1 syringe (50ml)

Size manikin size: approx. 18x8x10 (in)

Consumables and Replacement Parts

- 11393-010 face mask
- 11392-030 replacement tongue
- 11392-010 upper incisors

Recommended Sizes of Devices

Macintosh laryngoscope: blade size 3, 4 Tracheal tube: internal diameter

7.0mm, 7.5mm

Laryngeal Mask -Air-Q: size 3.5

-LMA Supreme: size4

-l-gel: size 4

Specifications are subject to change.