


New

Simulator Intravenous Arm IIa

IIa set is for USA and Canadian customers only

M50B-A



*The true-to-life IV training at budget friendly price.
Improved injection pad allows more repetitive training.
Practice at three veins: basilic, median antebrachial and cephalic veins.
Life-like vein palpation.*

Features Life-like vein palpation and realistic needle tip feeling.
Three veins for practice: Basilic, Median antebrachial and Cephalic veins.
Injection pads are easily replaceable.
No needle marks remain.
Simple and easy maintenance.

Specifications

Set includes:

- 1 arm model with supporter stand
- 1 injection pad for arm model
- 1 syringe
- 1 jar of coloring powder (red)
- 1 set of plastic beaker
- 1 spoon
- 1 injection bottle

Specifications are subject to change.

Replacement parts;

- 11267-210 replacement injection pad

