

# Patient Care Simulator "SAKURA"

M75



SAKURA is a multipurpose training model for teaching and learning hospital and community care procedures. The model can sit upright without support. Soft, lifelike materials and seamless limbs provide realistic hands-on training.

Life-like model with full articulation.

The model can sit upright without support.

Realistic articulating jaw and removable dentures facilitate oral care training.

Easy set up and simple clean up for catheterization training.

Catheterization fluids are retained in built-in reservoirs and are filled via ports in the hip.

Simulated stomach, held within and accessed via the chest wall, can be filled with fluids that can be extracted, or for checking nasogastric tube placement.

A cleaning system is incorporated into the simulator for flushing the urethra and rectum.

Removable basic intramuscular injection pads for arms, legs and hips.

Venipuncture trainer attachable to the arm is included.

Durable wig allowing hair combing / brushing, shampooing and blow drying.

Seamless arms and legs make the model more life-like.

An interchangeable male genitalia is available as an option (not included in the basic set).

Specifications are subject to change.

**For More Information Please Contact:**

**Becky Larson**

**Email: [Blarson@3-Dmed.com](mailto:Blarson@3-Dmed.com)**

**Direct: 952.447.6333**

**[www.3-Dmed.com](http://www.3-Dmed.com)**

## Patient Care Simulator "SAKURA"

### Training items

Body positioning  
 Patient handling  
 Removing and changing clothes  
 Bed bathing  
 Hair care  
 Oral & denture care  
 Tracheostomy care and management  
 Nasogastric tube replacement  
 Bag and mask replacement  
 Tracheal intubation care  
 Stomach pumping  
 Stoma care  
 Anal and rectal care  
 Intramuscular injection  
 Venipuncture and blood collection  
 Enema placement  
 Vaginal washing  
 Catheterization  
 Bedpan practice  
 Care of a deceased patient  
 Passive exercise

### Specifications

Set includes:  
 1 Sakura training model  
 fitted with a wig, dentures and a pair  
 of gluteal injection pads  
 1 pair of upper arm IM injection pads  
 1 pair of femoral IM injection pads  
 1 Strap-on Venipuncture Trainer  
 1 set of water bag and tube  
 2 kinds of drainage pump and tube ( 1 of each)  
 1 pair of pajamas  
 Specifications are subject to change.

Manikin size: 158Hcm, approx.18 kg  
 Packing size: 178 x 54 x 35 H cm, 27 kg

### Replacement and optional parts

11251-100 Male genitalia for M75  
 11251-010 Replacement wig for M75  
 11251-020 Replacement dentures for M75 (a set of upper & lower)  
 11251-030 Replacement a pair of shoulder injection pads for M75 (left & right)  
 11251-040 Replacement a pair of leg injection pads for M75 (left & right)  
 11251-050 Replacement a pair of hip injection pads for M75 (left & right)  
 11251-060 Replacement female genitalia for M75

