

Newborn Vital Signs Examination Simulator

M58



Available from December, 2007

Production Supervision
Keiko Takahashi, Professor,
College of Nursing Science, Fukui Prefecture University



Features

An audible, palpable and thermometric simulator replays vital signs of a newborn child, now get CE marked. Heart sounds, arterial pulsation at the wrists, neck and temples, and also rectal temperature can be checked. Heart sounds are recorded from actual patient and are synchronized with arterial pulses. All vital signs can be controlled.

Training items:

- Heart sounds
- Respiration
- Pulses: superficial temporal, carotid and radial arteries (bilateral)
- Body temperature (rectal thermometry)

Specifications

- Set Includes:
- 1 newborn manikin
 - 1 control box
 - 1 thermometer
 - 1 baby clothes
 - 1 carrying case
- Manikin Size: Approx. 50 cm
- Specifications are subject to change.



Palpation of pulses



Auscultation



Temperature taking

For more information contact:
Becky Larson
Email: Blarson@3-Dmed.com
Direct: 952.447.6333
www.3-Dmed.com

3-Dmed[®]
LEARNING THROUGH SIMULATION