## Strap-on Intramuscular Injection Trainer



The Strap-on Intramuscular Injection Trainer is designed for the practice between students and simulated patients. This trainer facilitates their understanding of the intramuscular injection process and the puncture feeling.

## Features:

- 1. Durable sensor is embedded to evaluate the procedure.
- 2. Both lamp and buzz show success or failur e in the procedure.
- 3. The water tank is located inside this trainer. It's compact and convenient.
- 4. This trainer is easier to drain without the waste water bag.
- 5. The skin and acr omion is similar to the feeling of a real human and it allows palpating the anatomical landmark.

## Strap-on Intramuscular Injection Trainer M67B

- 1 Base Flam
- Waste Water Tank
  Belt
- 1 Skin
  1 Soft Tissue Part
- 1 Soft Tissue Part 1 Muscular Part
- 1 Indication Box 1 AC adapter
- 1 Acromion Part

Specifications are subject to change.

Replacement Parts: 11297-210 Replacement muscular part with sensor 11297-220 Replacement skin 11297-230 Replacement soft tissue part



11297-210 Replacement muscular parts with sensor



11297-220 Replacement skin



11297-230 Replacement soft tissue parts

