Oral Care Model



"Effective oral care protects elderly people from pneumonia."

The importance of care has been paid attention to by many fields besides dental society, but in the medical, nursing and nursing care society as well.



Replacement of the Denture

Training Skills:

- Learn to observe oral cleanliness.
- Learn to maintain cleanliness of the denture and tongue.
- Denture replacement skills.



Life-like Sublingual



Maintenance for the Denture



Mirror Operation

Features:

- 1. This oral care model has a life-like textured sublingual and a replaceable denture.
- 2. This allows for practicing realistic oral cleaning procedure through use of simulated residue.
- **3.** Settings to accommodate differing patient postures necessary for practicing multiple procedures.
- **4.** The replaceable denture is realistic, allowing practical training of replacement and maintenance.



Oral Care Model M89

Set includes:

- 1 Oral Model with the Base Stand
- 1 Mirror
- 1 Tooth Brush
- 1 Denture Case
- 1 Simulated Residue

main body: Adult face model size: Life size material: Soft resin

Replacement Parts:

11277-020 Denture

Specifications are subject to change.