

Comprehensive training set with the platform of the maternal body for skills training in pelvic vaginal examination during labor and delivery assisstance.

Features

For all training skills

 Genital area made of innovative material that excels in elasticity, durability and resilience offers repetitive training opportunity.



Delivery assisstance

 Anatomically correct vaginal canal, inchial spine, pubic symphysis and coccyx that facilitate training with anatomical understanding.



 Elastic external genitalia, relistic shape of vaginal canal and the holding handle for the maternal body facilitate simulating a smooth delivery.

Cervical examination

 Interchangeable five stages of cervical dilation module can be installed with a single action.
Descent of the fetal head can also be set.



- Anatomically correct landmarks including the ischial spine and the public symphysis
- Anterior fontanelle to determine the fetal position.

Perineorrhaphy (perineal repair)

 Repetitive suture training for the first-degree perineal tear.



^{*} Suturing instruments are not included.

Cervical examination

Training Skills

- Cervical examination
 - ▶ Determination of station using landmarks of the ischial spine and the pubic joint to assess the descent of the body
 - ▶ Assessment of cervical dilation
- Speculum examination
- Procedures of observation and cervical cytology
- Procedures of urine catheterization



Pelvic vaginal examination



Procdures of urine catheterization

* For repracement parts, refer page 8.

Delivery assisstance

Training Skills

- Delivery assisstance for different delivery positions (dorsal, all-fours)
- Perineal protection
- Delivery of fetus (forceps delivery, vacuum extraction)

- Clamping, tying and cutting of the umbilical cord
- Delivery of placenta
- Insertion of a urine catheter



Perineal protection



Delivery assistance



Cutting the umbilical cord



Delivery of placenta

* For repracement parts, refer page 9.

Perineorrhaphy (perineal repair)

Training Skills

- Perineorrhaphy
 - ▶ Single interrupted suture
 - ▶ Vertical mattress suture



Cutting the umbilical cord

Specifications:

Common Compornents:	Parts for pelvic examination	Parts for Delivery Assistance	Parts for Perineorrhaphy
1 maternarl body	1 vagina unit	1 fetus	1 genitalia unit
1 genetalia unit (for pelvic examina-	1 set of cervical dilation module	1 placenta	(for perineorrhaphy)
tion and delivery assistance)	(5 variations)	3 umbilical cords for omphalotomy	(p
1 talcum powder	1 holder for the cervical dilation	1 vat	
3 lubricant (150ml)	module	5 velamen sheets	
1 instruction manual			