

M59/M63



# Newborn Bathing and Nursery Care Model



M59 Male

M63 Female

## Soft, seamless newborn manikin for skills training

Life-size 0-4-week-old infant, with the realistic weight and balance of newborn babies

### FEATURES

- 1 | Non-rigid neck flexibility
- 2 | Palpable anterior and posterior fontanel
- 3 | Bendable and stretchable limbs
- 4 | Seamless, one-piece waterproof structure
- 5 | Removable umbilical cord

### SKILLS

- | Bathing
- | Umbilical cord care
- | Bottle-feeding
- | Ear and nose care
- | Diaper and clothes changing
- | Physical measurement

### DESCRIPTIONS

#### SPECIFICATIONS

Size: 50 cm / 19.6 inch

Weight: 3 kg / 6.6 lbs

Head circumference: 33 cm / 13 inch

#### SET INCLUDES

- 1 | newborn manikin
- 1 | umbilical cord
- 1 | feeding bottle
- 1 | baby cloth
- 1 | instruction manual

#### MATERIALS

Soft resin  
Latex free

#### REPLACEMENT PARTS

11258-080 Umbilical cord

#### RECOMMENDED DEVICES

-

M58



# Newborn Vital Signs Examination Simulator

## Life-size newborn manikin for basic vital sign assessments

Product Supervision

Keiko Takahashi, professor  
College of Nursing Science,  
Fukui Prefecture University



### FEATURES

- 1 | Arterial pulse in wrists, neck and temples
- 2 | Heart sounds are recorded from actual patients and synchronized with arterial pulses
- 3 | All vital signs can be controlled with the control unit

### SKILLS

#### Vital signs examination

- | Heart sounds
- | Respiration
- | Pulses: superficial temporal, carotid and radial arteries (bilateral)
- | Body temperature (rectal thermometry)

### DESCRIPTIONS

#### SPECIFICATIONS

Size: 50cm / 19.6 inch

Weight: 1.1kg / 2.4 lbs

Power: AC100-240V

Power Consumption: 30W

#### SET INCLUDES

- 1 | newborn manikin
- 1 | control box
- 1 | thermometer
- 1 | baby cloth
- 1 | instruction manual

#### MATERIALS

Soft resin  
Latex free

#### REPLACEMENT PARTS

-

#### RECOMMENDED DEVICES

-