

MW2A



## Male Catheterization/Enema Simulator

MW52



## Female Catheterization/Enema Simulator II



MW2A



MW52

**High quality simulators for urethral catheterization, perineal care and enema practice**

Achieve better training results with our soft and realistic catheterization simulator that can also be used on SP or full body manikins.

**FEATURES**

- 1 | Soft and flexible genitalia allow proper cleaning and procedural preparation
- 2 | Successful insertion confirmation by fluid outflow
- 3 | Enema insertion in lateral position
- 4 | Detachable external genitalia

**SKILLS**

- |                                 |                                  |
|---------------------------------|----------------------------------|
| Intermittent catheterization    | Cleansing and aseptic technique  |
| Indwelling catheterization      | Catheter insertion and placement |
| Perineal care                   | Catheter Removal                 |
| Enema                           | Manual bladder compression       |
| Localization of urethral meatus |                                  |

Compact sets consist of a genitalia unit and accessories, for hybrid training and table top demonstration.

**MW2C** 11381-199-B

Male Genital Unit only

**MW2D** 11433-299-B

Female Genital Unit only



Enema Training



Pass of catheter can be observed through transparent bladder reservoir

**DESCRIPTIONS****SPECIFICATIONS**

Size: 34 x 33 x 24cm / 13.5 x 13 x 9.5 inch  
Weight: 3.5kg / 7.7lbs(M) 3.4kg / 7.5lbs(F)

**SET INCLUDES**

- 1 lower torso manikin
- 1 genital organ unit
- 3 valves (one is incorporated in the assembled unit)
- 1 irrigation bottle

**MATERIALS**

Soft resin  
Latex free

- 1 rectum tube
- 3 irrigation drainage tubes
- 1 lubricant
- 1 storage case
- 1 instruction manual

**REPLACEMENT PARTS**

- |  |                                       |
|--|---------------------------------------|
| 11381-100-01 clip (small) for urethral tube for MW2A | 11381-300 6 replacement valves        |
| 11381-100-02 clip (large) for MW2A/MW2B/MW2C         | 11433-010 1 replacement female skin   |
| 11433-000-01 clip (large) for MW52                   | 11433-020 2 replacement female meatus |
| 11381-100-03 urinary bladder reservoir               | 11381-800 lubricant                   |
| 11381-100-04 rectum unit                             | 11381-410 1 replacement male skin     |

- |                                     |
|-------------------------------------|
| 11381-510 2 replacement male meatus |
|-------------------------------------|

**RECOMMENDED DEVICES**

\* Catheters are not included.