



**External genitalia made of the new material and anatomically correct vaginal canal provides true-to-life hands-on training in cervical examination to determine Bishop score from the final stage of pregnancy to the first stage of the labor.**

### Features

1. Genital area made of innovative material that excel in elasticity, durability and resilience offers repetitive training opportunity.
2. Interchangeable five stages of cervical dilation module can be installed with a single action. Descent of the fetal head can also be set.
3. Anatomically correct landmarks including the ischial spine and the pubic symphysis
4. Posterior fontanelle to determine the fetal position

### Training Skills


- Cervical examination
  - ▶ Determination of bishop score (cervix score) using landmarks of the ischial spines and the pubic symphysis
  - ▶ Assessment of cervical dilation
  - ▶ Procedures of urine catheterization
  - ▶ Procedures cervical cytology
- Speculum examination

### Key Features

**Eternal genitalia that excels in elasticity, resilience and durability**



**Ischial spines and the pubic symphysis to determine fetal position**



Station zero, the line that connect ischial spines can be recognized.



Cervical dilation module can be attached by a single action.



Descent of the fetus can be set by the lever on the holder. (from -3 to +6)

### Specifications :

Set Includes:		Manikin Size:	Replacement Parts:
1 maternal body	1 holder for the cervical dilation module	W74 x D44 x H33.5 cm	(11415-010)
1 genitalia unit	3 lubricant (150 ml)		1 set of cervical dilation module (5 variations)
1 cervix unit	1 talcum powder		(11415-020) 3 lubricant (150 ml)
1 set of cervicildilation module (5 variations)	1 instruction manual		(11415-030)
			1 genitalia unit (for pelvic examination and delivery assistance)