

Ultrasound Lumbar Puncture Simulator

This Lumbar Puncture/Epidural Simulator is an ultrasound compatible simulator that includes the lumbar vertebrae, iliac crest, spinous process, superior articular process, transverse process, epidural space and dura mater. In addition, the main body is common with Lumbar Puncture Simulator and Fluoroscopy Lumbar Phantom. Therefore, all lumbar puncture replacement blocks can be used with this model.

Point 01 Multi device Trainer

Replacement parts are durable for multiple procedures.



Ultrasound Guided Lumbar Puncture




Fluoroscopy Guided Lumbar Puncture



Landmark Guided Lumbar Puncture

Point 02 Positioning

This simulator can be positioned in the upright or lateral position.



- Transparent block makes the needle path visible.
- Anatomically correct body and block provide realistic ultrasound image, and platform for landmark procedure.
- Latex free

Training Skills

1



- 1 Ultrasound guided lumbar puncture
2 Ultrasound guided epidural procedure

2



3



- 3 Landmark guided lumbar puncture
4 Landmark guided epidural procedure

4

5



- 5 CSF collection and CSF pressure measurement