# CVC Insertion Simulator III

## **CVC** III facilitates a sequence of uninterrupted simulation experience for specialist training

From aseptic field preparation to final catheter placement









### **FEATURES**

- 1 | Insertion practice is repeatable to obtain the best technique for patient safety.
- 2 | Both landmark guided insertion and ultrasound guided techniques can be practiced.
- 3 | The frictionless tissue-like simulation pad to facilitate a realistic sensation of the CVC placement.

### SKILLS

Take complete steps to a successful CVC insertion


| Perform hand hygiene practices and prepare aseptic field

| Decide on insertion site via ultrasound guidance or observation and palpation of landmarks

| Insertion under ultrasound guidance or via landmark technique | Confirm flashback

| Advance guidewire and sheath

Insert catheter





Anatomically accurate features including 3 insertion sites and rib cage

## **DESCRIPTIONS**

SPECIFICATIONS Size: W40 x D20 x H 34 cm W15.7 x D8 x H13.5 inch

#### SET INCLUDES

- 1 male upper torso manikin
- CVC placement pads
- 1 skin for cannulation training transparent anatomical block
- 2 artery tube (transparent): 2 type
- 2 vein tubes (blue): 2 type
- vein pipes
- pulsation unit (for landmark technique)
- irrigation bottle (for UG technique)

MATERIALS Soft resin, hydrogel, hard resin Latex free

- 1 sample needle
- syringes (2.5ml, 50ml)
- 1 introductory ultrasound training block REAL VESSELS

REPLACEMENT PARTS

11347-210 2 introductory ultrasound training blocks REAL VESSEL

11347-240 2 CVC placement pads 11347-260 1 pulsation unit

11347-320 irrigation bottle

Ultrasound scanner and needles are not included.